

Porter Instrument's Nitrous Nasal Hood Solutions Easily Retrofit To Any Manufacturer's Breathing Circuit

PORTER
The Trusted Name In Nitrous

PORTER AUTOCLAVABLE DOUBLE MASK



- The only autoclavable double mask in the industry.
- The soft and pliable inner liner creates a tight, yet comfortable seal around the patient's nose, assuring superior scavenging efficiency.
- Available in both Adult and Pediatric sizes.
- Easily retrofits Accutron™ or Belmed breathing circuits, no adapter needed.



Porter Double Mask
Autoclavable

DYNAMITE SINGLE-USE NASAL HOODS



- Single-use disposable nasal hood with dinosaur sticker.
- Soft, pliable material allows for a comfortable fit.
- 3 size options and 5 scents to choose from.
- Easily retrofits to Accutron™ and Belmed breathing circuits with the Universal Conversion Kit.




Dynomite Single-Use Scented
Nasal Hoods



Universal Conversion Kit

SILHOUETTE LOW PROFILE NASAL HOOD



- Single-use nasal hood and breathing circuit **ALL IN ONE!**
- Disposable mask and circuit for improved infection control.
- Unobstructed access to the oral cavity.
- Four sizes for a perfect fit on all patients.
- Easily retrofits to all flowmeters with Silhouette Retro Package.



Silhouette Single-Use Nasal Hood
And Breathing Circuit



Silhouette
Retro Kit

PORTER AUTOCLAVABLE NASAL HOODS



5054C

Porter Adult Double Mask w/ 1 Inner Mask Liner

\$351

5054D

Porter Pediatric Double Mask Hood Kit w/ 1 Inner Mask Liner

\$351

Retrofitting to Porter is easy!



Accutron Clearview™



1. Disconnect old mask from breathing circuit.



2. Connect new Porter mask to breathing circuit.



Accutron PIP+™



1. Disconnect old mask from breathing circuit.



2. Connect new Porter mask to breathing circuit.


DYNAMITE 24 PACK BAG

| | | |
|----------|-------------------------|-------|
| 91316461 | Bubble Gum Small | \$146 |
| 91316468 | Bubble Gum Medium | \$146 |
| 91316475 | Bubble Gum Large | \$146 |
| 91316462 | Strawberry Small | \$146 |
| 91316469 | Strawberry Medium | \$146 |
| 91316476 | Strawberry Large | \$146 |
| 91316463 | Orange Small | \$146 |
| 91316470 | Orange Medium | \$146 |
| 91316477 | Orange Large | \$146 |
| 91316464 | Vanilla Small | \$146 |
| 91316471 | Vanilla Medium | \$146 |
| 91316478 | Vanilla Large | \$146 |
| 91316465 | Plain Small | \$146 |
| 91316472 | Plain Medium | \$146 |
| 91316479 | Plain Large | \$146 |
| 91316466 | Assorted Scents Small | \$146 |
| 91316473 | Assorted Scents Medium | \$146 |
| 91316480 | Assorted Scents Large | \$146 |
| 91316481 | Assorted Scents Various | \$146 |

DynoMite



Retrofitting to Dynomite is easy!



Accutron Clearview™



1. Disconnect old mask from breathing circuit.



2. Connect Universal Conversion Kit to breathing circuit.



Accutron PIP+™



1. Disconnect old mask from breathing circuit.



2. Connect Universal Conversion Kit to breathing circuit.



UNIVERSAL CONVERSION KIT

91515142 DynoMite Nasal Hood Universal Conversion Kit \$114

SILHOUETTE ALL IN ONE SINGLE-USE NASAL HOODS/BREATHING CIRCUITS

| | | |
|---------------------|--|--------------|
| SIL2-PED0-12 | Pediatric Silhouette – 12 Pack | \$114 |
| SIL2-SM-12 | Small Silhouette – 12 Pack | \$114 |
| SIL2-MED-12 | Medium Silhouette – 12 Pack | \$124 |
| SIL2-LG-12 | Large Silhouette – 12 Pack | \$124 |
| SIL2-PED0-24 | Pediatric Silhouette – 24 Pack | \$213 |
| SIL2-SM-24 | Small Silhouette – 24 Pack | \$213 |
| SIL2-MED-24 | Medium Silhouette – 24 Pack | \$239 |
| SIL2-LG-24 | Large Silhouette – 24 Pack | \$239 |
| SIL2-VAR-4x3 | Variety 12 Pack: 3 of each size - Pediatric, Sml, Med, Large | \$124 |

Porter Silhouette can be retrofitted to any flow meter on the market. [Here's all you need!](#)



SIL2-RETRO-PKG

Flowmeter Connection Kit (order one for each flowmeter)

\$99

SIL2-VAR-4x3

Variety 12 Pack

\$124



Click the link below for more information:

www.PorterInstrument.com/NasalHoods

Or contact us directly at:

porternitrous@parker.com

215-723-4000

PORTER
The Trusted Name In Nitrous



Accutron™ is a registered trademark of Accutron, Inc., a division of Crosstex, a Cantel Medical Company.